



Simplify data transfer from scales and balances with Adam Connect wireless connectivity. The plug-and-play Adam Connect Bluetooth module works with compatible weighing devices to transmit data up to 35 ft. away to PCs or printers using the Adam Connect USB PC Connector (for data collection using the AdamDU software or other LIMS systems) and Printer Connector modules. For a customized setup, select only the receiver modules you need. The Bluetooth module is required with each configuration.

Plug-and-play wireless connectivity for data communication between scales and balances with PCs and printers



Item Number.....	Description	Item Code
8000011571	Bluetooth Sender	BTM-01
8000011587	Bluetooth-to-Printer Receiver	BTM-02
3074010507	Bluetooth-to-PC Receiver	BTM-03

About ADAM:

For more than 50 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Adam Equipment follows a policy of continuous improvement and reserves the right to change specification and pricing without notice. Adam Equipment's standard terms and conditions of sale apply. Errors and omissions excepted.