



Simplify data transfer from scales and balances with Adam Connect wireless connectivity. The plug-and-play Adam Connect Bluetooth module works with compatible weighing devices to transmit data up to 35 ft. away to PCs or printers using the Adam Connect USB PC Connector (for data collection using the AdamDU software or other LIMS systems) and Printer Connector modules. For a customized setup, select only the receiver modules you need. The Bluetooth module is required with each configuration.

Plug-and-play wireless connectivity for data communication between scales and balances with PCs and printers



Item Number	Description	Item Code
8000011571	Bluetooth Sender	BTM-01
8000011587	Bluetooth-to-Printer Receiver	BTM-02
3074010507	Bluetooth-to-PC Receiver	BTM-03

About ADAM:

For more than 50 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Adam Equipment follows a policy of continuous improvement and reserves the right to change specification and pricing without notice. Adam Equipment’s standard terms and conditions of sale apply. Errors and omissions excepted.

Product Details

Model	AE 505
Capacity	Selectable
Readability	Selectable
Weighing Units	g, kg, lb, lb:oz, oz, ct, dr, mm, tL,T, T
Interface	RS-232, USB, Ethernet RJ45, Wi-Fi, Bluetooth
Calibration	External
Display	Touchscreen with #mm-high digits
Power Supply	12VDC 3A adapter, Internal rechargeable battery
Operating Temperature	32° to 104°F / 0° to 40°C
Housing	ABS Plastic
Overall Dimensions (WxDxH)	11.8" x 4.9" x 6.7" / 300 x 170 x 125mm

Accessories

Item NumberDescription

400005615..... Labels — 800 per roll
 3014011014..... RS-232 cable 9M-9F
 3074010507..... RS-232 to USB adapter
 600002028..... AdamDU Data Collection Program
 700660290..... Calibration certificate

Key Accessories



Labels

A roll of 800 sticky labels intended for use with our AE 503, 504 and 505 Label Printing Indicators.



AdamDU Data Collection Program

Adam DU quickly collects information from up to eight Adam balances simultaneously. An easy-to-use interface with a main command window allows the user to view and manage data capture sessions.